Spin Noise Spectroscopy of Rubidium

David Blasing

Professor: Dr. Yong P. Chen Advisor: Sourav Dutta

- Review of Pertinent Quantum Mechanics
- Hyperfine Structure of Rubidium 85
- Faraday Rotation
- Spin Noise
- Future Goals

Review of Pertinent Quantum Mechanics

- It is possible to excite an atom (or system) to a higher energy level using a tuned laser (E=hv)
- The system then transitions to a lower energy state by emitting a photon with energy equal to the energy difference between states
 - E_{photon} = E_{excited} E_{ground}
 - The decay rate, γ, of rubidium is about 3.6x10[^]7/s transitions per second. Alternatively, the average lifetime of the excited states studied here will be about 28 ns ¹
- Emission is spontaneous (random direction)

Transition rate is a ballpark for the rubidium levels reached by the 780nm laser


Review of Pertinent Quantum Mechanics

- ${\mbox{\ensuremath{\bullet}}} F$ is the magnetic quantum number, and m_f is the projection of J onto the quantization axis
- •The presence of an external magnetic field breaks the degeneracy for the different values of m_f. This arises from the the potential energy -µ•B.
- •The weak field zeeman effect predicts these energies as $\Delta E = g_f \mu_B m_f B$
- •Right handed circularly polarized light σ_+ carriers momentum +h_{bar} causing transitions from m_f to m_f+1 while σ_- (left handed) carries momentum $-h_{bar}$, causing transitions from m_f to m_f-1


"If you are not completely confused by quantum mechanics, you do not understand it." -Theoretical Physicist John Wheeler


 G_j is the lande G factor, μ_B is the bohr magneton, m_j is the z-component of the total angular momentum, and B is the applied magnetic field. The magnetic field associated with the fine structure correction is ~3000 Gauss.


- Review of Pertinent Quantum Mechanics
- Hyperfine Structure of Rubidium 85
- Faraday Rotation
- Spin Noise
- Future Goals



Fine structure results from spin-orbit coupling of the electron, while hyperfine structure results from spin of the nucleus interacting with the spin of the electron.









Fine structure results from spin-orbit coupling of the electron, while hyperfine structure results from spin of the nucleus interacting with the spin of the electron.

- Review of Pertinent Quantum Mechanics
- Hyperfine Structure of Rubidium 85
- Faraday Rotation
- Spin Noise
- Future Goals



Positive faraday constant means counter clockwise rotation when then direction of the magnetic field and laser are parallel, and clockwise when they are anti parallel.




Positive faraday constant means counter clockwise rotation when then direction of the magnetic field and laser are parallel, and clockwise when they are anti parallel.


In the linearly polarized light equation, the k's are different because $\boldsymbol{\lambda}$ is frequency dependent


Faraday Rotation

- •The index of refraction is depends on the frequency of light •n(v)-1=n_q \sigma_0 (\gamma_0/4\pi)^2/((v_0-v)^2+(\gamma_0/4\pi)^2) ^2
- •The Zeeman Effect
 - •A weak magnetic field associates each m_f sublevels with a different energy
 - •Right and left circularly polarized light cause different transitions within these sublevels
 - •Right/left circular have different resonant energies, and therefore different indices of refraction


 n_g is the number density of atoms σ_0 is the integrated cross section of interaction $\gamma_0/2\pi$ is the natural line width of the transition ν_0 is the resonant frequency









The peaks are both from Rb 85, f=2,3 to excited states



Larger Peak is Rb 85 while the smaller peak is Rb 87

- Review of Pertinent Quantum Mechanics
- Hyperfine Structure of Rubidium 85
- Faraday Rotation
- Spin Noise
- Future Goals





$$\langle \delta \theta_{\rm F}^2 \rangle^{1/2} = \frac{e^2 f \beta}{4 m c \Delta} \sqrt{\frac{N_0 L}{A}}$$

Variable Parameters:


- •Homogeneity of the magnetic field
- •∆ (detuning from resonance)
- •N₀ (density of atoms)
- •L (interaction length)
- •A (cross sectional area of laser beam)


$$\langle \delta \theta_{\rm F}^2 \rangle^{1/2} = \frac{e^2 f \beta}{4 m c \Delta} \sqrt{\frac{N_0 L}{A}}$$

- •Magnetic field is 1.85 Gauss with no more than %10 variation over Rubidium cell
- •Optimized the detuning of the laser based on previous results, $\Delta{\sim}20\mbox{Ghz}$ from the \mbox{D}_2 transition
- •10 cm length of Rubidium cell is longer than previous experiments
- •Heat tape increases density of rubidium to ~10⁷ atoms/mm³
- •Circular aperture and lenses shrink beam area ~.12 mm²

- •Spin noise should be occurring at the larmor frequency, $g_f \mu_b B/h$
- •Larmor frequency for ⁸⁵Rb~.87Mhz, ⁸⁷Rb~1.32Mhz
- •Currently, we see now difference between on the spectrum analyzer when the magnetic field is on vs. magnetic field off
- •The last portion of the R.E.U. was dedicated to maximizing the parameters that govern spin noise, and trouble shooting the amplification circuit

H=planks constant, μ_b =planks constant, B is magnetic field, g_f =ground state g-factor.






 R_t determines the bandwidth and risetime of the photodiodes. $R_t \! = \! 1k\Omega \! \sim \! 4Mhz.$

The gain, R_f/R_i , determines the bandwidth of the gain stage. For OPA228, R_f/R_i =3 means bandwidth~2.6 Mhz.

Total transimpedance gain ~ 3V/mA.

Previous groups had transimpedance gains of ~40V/mA



With the 3 gain stage circuit, we attempted to find spin noise one last time

- In the spectrum analyzer around 1Mhz
 •Noise floor was ~500μV with fluctuations of 250μV.
 •Two peaks of 3.1mV and 7.7mV with fluctuations of 1.5mV and 2.5mV
- •Spin noise will be in the nanovolt region

- Review of Pertinent Quantum Mechanics
- Hyperfine Structure of Rubidium 85
- Faraday Rotation
- Spin Noise
- Future Goals

Future Goals

Two main areas for future work

- 1. Refining the controllable parameters for spin noise
 - Creating a more static magnetic field, or using a smaller rubidium cell
 - More elaborate set of optics to shrink the beam area even further
 - Cell oven to increase max temperature
- 2. Engineering a low noise, high bandwidth, high gain circuit
 - Current is on the scale of picoamps or lower, and is oscillating around 1mhz
 - Commercially balanced, low noise detectors are available, but are pricey
 - Examining each circuit element in the spectrum analyzer to see its noise contribution, and selecting lower noise alternatives if neccesary

Acknowledgements

I would like to thank Dr. Chen, my adviser Sourav Dutta, the grad students the lab and the physics department for allowing, and helping me contribute to A.M.O. research done at Purdue.

Works Cited

- 1) University of Florida. Saturated Absorption Spectroscopy.
 Advanced Laboratory Manual.
- 2) Waitukaitis, Scott. Resonant Faraday Rotation In A Hot Lithium Vapor. University of Arizona.
- 3) Crooker et. al. Spectroscopy of spontaneous spin noise as a probe of spin dynamics and magnetic resonance, *Nature 341*, 49-52 (2 September 2004).
- 4) Mills, A.. Nonlinear ground-state pump-probe spectroscopy in ultracold rubidium: Raman-coupled dressed state spectroscopy. *Diss. Purdue University 2007*. Dissertations & Theses @ CIC Institutions, ProQuest. Web. 26 Jul. 2010