# Biophysics: Diffusion through Biological Tissue

Andrew McElroy
Summer REU 2011
Purdue University


## Diffusion

- What is it?
  - The net movement of particles from an area of high concentration to an area of low concentration.
  - Constant random motion:
    - Brownian motion
    - "Random Walk"



# Why study it?

- Diffusion is well defined for simple systems.
- There is a lot to be learned about more complex systems.














## Research Applications

- Understanding of the transport properties of a cell
  - Perception of a cell:
    - Water sack with stuff floating around inside
    - Not true
    - By studying diffusion we can gain a better knowledge of the true characteristics of cells

# Research Applications

- Cancer treatment
  - Nanoparticles are becoming a hot bed of research.
    - Supposed better drug transportation and targeting.




# Experiment

- Cell Design needs:
  - Two connected chambers.
    - High concentration and low concentration
  - Ability to hold a sample between the chambers
  - Interchangeability
  - Way to measure the concentration changes.

#### Cell:

Two connected chambers





- Holding a sample/interchangeability
  - Compress a plug of a sample between the cells




# **Making Measurements**

- Needed a way to measure concentration changes without disturbing the system.
  - Use fluorescing dyes




 Rhodamine and other fluorophores can easily be detected and can be used to tag other larger molecules

### • System works on our nifty scanner:











#### Problem

- After all the simple problems (leaking, etc)
  appeared to be solved in the system, the
  experiments were still not acing like they
  should.
- There was consistently a bulk amount of dye entering the system early in the experiment.




#### Problem

- Inherent problem in diffusion measurements:
  - Diffusion is a very slow process that is often drowned out by convection currents and other forces







# **Problem Solving**



-Let it sit and added membranes to both ends of the "tunnel" to hopefully stop these convection currents











# **Problem Solving**



Need to change our approach

# System 2.0






#### Density issues?

- Solved by making the bottom solution more dense than the dye



## DATA!




Rate of change (slope):

trial 1: .43uM/hr

trial 2: .41uM/hr

# Steady State Model



$$\frac{\partial C}{\partial t} = \frac{DAC_0}{VL} \left( \frac{1}{1 + \frac{2\Delta x}{Lf}} \right)$$

 This gives an expected rate value which makes the measured value factor of 3 too slow.... Not sure why

#### Conclusion

 Successfully made a cell which can make diffusion measurements through biological tissue and gells

Begin measuring the diffusion rate through agarose gel

 Move on to measuring diffusion through biological tissue